

Summary of Solutions Lab Calculations



Calculations *Show all your work*

You will need the molar mass of NaCl

1. Calculate the volume of 2.000 M NaCl solution needed to obtain 1.000 gram of NaCl.

$2\text{ M} = 2\text{ moles solute per liter Or } 2\text{ moles solute per } 1000\text{ mL}$

$$1.000\text{ g NaCl} \times \frac{1\text{ mole NaCl}}{\text{Molar Mass g NaCl}} \times \frac{1000\text{ mL}}{2.000\text{ moles}} = ?\text{ mL of } 2.000\text{ M NaCl solution}$$

2. Calculate (by subtraction) the starting mass of NaCl solution (weighing by difference):

Mass of evaporating dish, watch glass, and 2.000 M NaCl solution

- Mass of evaporating dish, watch glass
Mass of NaCl solution

3. Calculate the mass of isolated solid NaCl (weighing by difference):

Mass of evaporating dish, watch glass, and white solid in the evaporating dish

- Mass of evaporating dish, watch glass
Mass of NaCl (the actual or experimental yield)

4. Calculate the % yield for the experiment

$$\% \text{ Yield} = \frac{\text{Mass isolated NaCl (from 3 above)}}{\text{Mass theoretical NaCl (1.00 g)}} \times 100$$

5. Calculate the % error for the experiment

$$\% \text{ Error} = \frac{\text{Actual yield (g)} - \text{Theoretical yield (g)}}{\text{Theoretical yield (g)}} \times 100$$

