



Acids $\xleftrightarrow{\text{destroy}}$ Bases


acid $\text{HCl} + \text{NaOH}$ base

$\text{NaCl} + \text{H}_2\text{O}$

salt water

Titrations



Copyright Larry P. Taylor, Ph.D. All Rights Reserved LPT

Titration Terms

Titration
Controlled addition of a liquid into a vessel to measure the volume that reacts with a substance already in the vessel


Indicators
substances that change color to signal when to stop a titration
Organic dyes whose color is sensitive to pH

Endpoint
point in a titration when the indicator changes color

Standard Solution
solution of known concentration used in a titration

Neutralization
double replacement reaction: an acid and a base react to form water and a salt

$\text{Acid} + \text{Base} \rightarrow \text{Salt} + \text{Water}$





Copyright Larry P. Taylor, Ph.D. All Rights Reserved LPT

Titration Experiment

Indicator Color	0	1	2	3	4	5	6	7	8	9	10	11	12	13
Universal Indicator	Red	Orange	Yellow	Green	Blue	Blue	Blue	Blue	Blue	Blue	Blue	Blue	Blue	Blue
Phenolphthalein	Colorless	Colorless	Colorless	Colorless	Colorless	Colorless	Colorless	Colorless	Colorless	Colorless	Colorless	Colorless	Colorless	Colorless
Methyl Orange	Red	Red	Red	Red	Red	Red	Red	Red	Red	Red	Red	Red	Red	Red
Bromothymol Blue	Red	Orange	Yellow	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green
Phenol Red	Red	Red	Red	Red	Red	Red	Red	Red	Red	Red	Red	Red	Red	Red
Thymol Blue	Red	Orange	Yellow	Green	Blue	Blue	Blue	Blue	Blue	Blue	Blue	Blue	Blue	Blue



Flask contains 0.1M strong base.

Buret contains 0.1M strong acid.

Flask contains 0.1M strong acid.

Buret contains 0.1M strong base.

End Point

Copyright Larry P. Taylor, Ph.D. All Rights Reserved LPT

Titration of Strong vs. Weak Acids

Strong Acid (HCl)

Titration of strong acid (25 cm³) with strong base.

Weak Acid (HOAc)

Titration of weak acid (25 cm³) with strong base.

Same End-point

Copyright Larry P. Taylor, Ph.D. All Rights Reserved LPT

Phenolphthalein

phenolphthalein
pH < 7.2

phenolphthalein
pH > 8.2

One of the most common indicators used

Laxative

C.S.I. = preliminary test for blood

Kastle-Meyer Spot Test

Phenolphthalein plus sample

Add H₂O₂

Hemoglobin present oxidizes to pink form

OH⁻ attacks acid form → changes structure

Acid form: colorless

Basic form: magenta

Copyright Larry P. Taylor, Ph.D. All Rights Reserved LPT

Titration Experiment

Key to "Titrations"

At end point:

Moles standard added = moles unknown present

Copyright Larry P. Taylor, Ph.D. All Rights Reserved LPT

Neutralization Reactions: Solution Stoichiometry

At Endpoint: moles added = moles unknown



All titration problems solved the same way:

Balance the chemical Reaction

Determine moles present in standard solution (moles / L x L)

Use reaction coefficients ("per expression") to get moles unknown

Convert moles of unknown

to solution concentration (molarity)

to grams present

to gas volume




Copyright Larry P. Taylor, Ph.D. All Rights Reserved

LPT

Titration of Acetic Acid with Sodium Hydroxide

Vinegar = dilute solution of acetic acid (CH_3COOH or $\text{HC}_2\text{H}_3\text{O}_2$ or HOAc) in water

The acetic acid will react with a base such as sodium hydroxide (NaOH)



Acid + Base \rightarrow Salt + Water


Problems

Typically prepared NaOH solution is not well characterized:

Solid NaOH readily absorbs moisture from the air.

Initial weighing error

Atmospheric CO_2 reacts with water to make carbonic acid.



Acid reacts with some NaOH \rightarrow lowers concentration of the NaOH



Copyright Larry P. Taylor, Ph.D. All Rights Reserved


LPT

Solution



Concentration of the NaOH is determined by titration with known strong acid

Titration NaOH with standard acid solution \rightarrow "standardization"



Standardization gives accurate value of NaOH solution.

Standardized NaOH titrated against unknown acid

At Endpoint: Moles standard NaOH = moles unknown

Standard NaOH then used to titrate acid concentrations



Copyright Larry P. Taylor, Ph.D. All Rights Reserved


LPT

Keep Your End-Point Indicator Barely Visible



Copyright Larry P. Taylor, Ph.D. All Rights Reserved

LPT



Titration Lab



Copyright Larry P. Taylor, Ph.D. All Rights Reserved

LPT

Online Lab



Purpose:
To determine the molarity of acetic acid


Since NaOH concentration can vary with time, the lab is done in 2 parts:

1. Determine the molarity of NaOH solution using standard HCl
2. Using the standardized NaOH solution, determine the molarity of HOAc

For Measuring Acid


For Dispensing Base

Reaction Happens Here



Copyright Larry P. Taylor, Ph.D. All Rights Reserved

LPT



Data

Titration Volumes

Calculations

First Titration: Molarity NaOH

Molarity of HCL (from instructor): 0.7664

$$\text{HCl (mL)} \times \frac{\text{HCl Moles}}{1000 \text{ mL}} \times \frac{1 \text{ mole NaOH}}{1 \text{ mole HCl}} \times \frac{1}{\text{NaOH added from buret (L)}} = \text{Molarity NaOH}$$

$$25.00 \text{ HCl mL} \times \frac{0.7664}{1000 \text{ mL}} \times \frac{1 \text{ mole NaOH}}{1 \text{ mole HCl}} \times \frac{1}{0.02398 \text{ L}} = \text{Molarity NaOH}$$

Second Titration: Molarity HOAc

Molarity NaOH (from Titration 1)

$$\text{NaOH added (mL)} \times \frac{\text{NaOH Moles}}{1000 \text{ mL}} \times \frac{1 \text{ mole HOAc}}{1 \text{ mole NaOH}} \times \frac{1}{\text{HOAc in Erlenmeyer (L)}} = \text{Molarity HOAc}$$

$$28.48 \text{ mL} \times \frac{0.7990 \text{ moles}}{1000 \text{ mL}} \times \frac{1 \text{ mole HOAc}}{1 \text{ mole NaOH}} \times \frac{1}{0.02500 \text{ L}} = \text{Molarity HOAc}$$

Copyright Larry P. Taylor, Ph.D. All Rights Reserved

Results

Fill In table

Conclusion

Fill-In molarity value

Nothing needed but that value

Let's Boldly Go Explore Today's Lab



Copyright Larry P. Taylor, Ph.D. All Rights Reserved

LPT
